

Supplementary Material

Slide-to-Slide Tissue Transfer and Array Assembly from Limited Samples for Comprehensive Molecular Profiling

Stephanie E. Weissinger^{1,2}, N. Zeke Georgantas³, Julia C. Thierauf³, Rebecca Pellerin³, Emma Gardecki³, Stephanie Kühlinger¹, Lauren L. Ritterhouse³, Peter Möller², and Jochen K. Lennerz³

¹*Institute of Pathology, Alb Fils Clinics GmbH, Göppingen, Germany*


²*Institute of Pathology, University Hospital Ulm, Ulm, Germany*

³*Center for Integrated Diagnostics, Department of Pathology, Massachusetts General Hospital/Harvard Medical School, Boston, USA*


Content

Supplementary Figure 1: Configuration of the STS components.....	2
Supplementary Figure 2: Examples of STS Transfer in the context of incident management.....	3
Supplementary Figure 3: Competency Assessment and Technical Success Rate of the STS Array Assembly Technique.....	4
Supplementary Video Description.....	5

Supplementary Figure 1




Supplementary Figure 2



Supplementary Figure 2. Example of STS Transfer in the context of incident management. (a) The suspected contaminant on the H&E slide is present in multiple levels; contaminant was not present in deeper levels (not shown). (b) After lifting and transfer to a charged slide, additional immunohistochemical workup clarified distinct origin. *Abbreviations:* H&E, haematoxylin & eosin; IHC, immunohistochemistry; STS, slide-to-slide.

Supplementary Figure 3



Supplementary Figure 3. Competency Assessment and Technical Success Rate of the STS Array Assembly Technique. To practice the STS technique or to assess manual competency in performing the lift–split–mount steps (steps 4/5/6), we use a replaceable H&E section. (a) After coverslip lifting and xylene immersion the slide is covered with methacrylate. After drying, the slide is immersed in water overnight (not shown). Subsequent steps include lifting, splitting and mounting (see [Supplementary Video](#)). (b) The assembled STS array allows assessment of the manual competency in mounting the individual pieces ('tiles'; manual completion rate). (c) The same array can be used to assess technical success rate; shown here as an immunohistochemical stain for cytokeratin 7. The technical success rate is affected by the principal antigen retrieval technique (and sub-conditions) as well as the counterstaining and final coverslipping methods. The overall performance of the STS array technique is a combination of manual competency during assembly and technical/analytical success. For educational purposes we depict an STS array with 2 complete tissue dropouts, partial losses, and folds. *Abbreviations:* H&E, haematoxylin & eosin; IHC, immunohistochemistry; STS, slide-to-slide transfer

Supplementary Video Description

Clip that outlines all steps of the slide-to-slide (STS) transfer technique including a practice exercise for competency assessment before working with precious tissues.

Length: 2:56 min

Format: .mp4